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## SAMPLING PROTOCOL FOR HIV VIRAL LOAD

### **Blood sampling:**

9 ml (for adult) / 3,5 ml (for child) EDTA-blood – mix well. Use plastic blood tubes for the sampling (**glass tubes will be refused**).

Store blood at room temperature (20-25°C).

### **Request:**

There must be a completed and signed request form for each sample.

On the sample: give date and hour of the blood sampling.

### **Transport blood:**

Via a courier (not by normal mailing).

Transport immediately to ITM (**within 24 hours after blood sampling**).

Delivery on weekdays from 8:30am-3:00pm.

### **Plasma isolation in your own laboratory:**

If the blood cannot be delivered within 24 hours, plasma can be isolated in your own laboratory, in accordance with the following procedure.

1. Centrifuge the blood tube during 20 minutes at 1100xg (approx. 2200-2500rpm) at room temperature.
2. Divide the plasma (3-5ml) in about 2-3 micro tubes of 1.5ml-2ml with a minimum 1.2ml plasma per tube (we need at least 2,5ml plasma for the analysis).

Use micro tubes, resistant to -80°C -type Sarstedt 72.694.005 or type Greiner 722361.

**Larger tubes are denied.**

3. Freeze the micro tubes at -20°C (± 5°C) or -80°C (±15°C), **in upright position. The plasma should be fully frozen before transport.**

### **Transport of frozen plasma:**

Via courier.

Transport on dry ice or with cool packs.

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### **Used test and test limits:**

The Hiv-1 viral load is measured with the Cobas 4800 HIV-1 test, with a test range of 20 – 10,000,000 RNA copies/ml.

The Hiv-2 viral load is outsourced to another ARL. They use an in house method with a test range of 10 – 5,000,000 RNA copies/ml.

All tests happen under BELAC accreditation.

### **How soon can you expect results:**

- HIV-1 Viral load: ca. 3 weekdays
- HIV-2 Viral load: ca. 15 weekdays

### **For more information:**

Director ARL – Dorien Van den Bossche – Phone: 03/345 55 97 – E-mail: [dvandenbossche@itg.be](mailto:dvandenbossche@itg.be)

Study- and Reference coordinator ARL – Fien Vanroye – Phone: 03/247 66 14 – E-mail: [fvanroye@itg.be](mailto:fvanroye@itg.be)